Tel-O-Graph® GSM

Tel-O-Graph® GSM plus

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DE

Operating instructions

EN

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ES

Mode d'emploi

FR

Tel-O-Graph® GSM

Tel-O-Graph® GSM plus

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DE

Tel-O-Graph® GSM
Tel-O-Graph® GSM plus

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The contents are not permitted to be duplicated or published without written authorization by I.E.M. GmbH.

Revision: EN_1.0 / 21.04.2016

© I.E.M. GmbH 2016. All rights reserved.

Content

1	Introduction	46	4.3	Switching the monitor on/off	57
2	Description of device	47	5	Measuring blood pressure and pulse	58
2.1	The monitor	47	5.1	Before measuring	58
2.2	Cuff	48	5.2	Armmanschette anlegen	58
2.3	Display	49	5.3	Correct posture	60
2.4	Content	50	5.4	Measuring	61
2.5	Clinical Validation	50	5.5	Stopping the measurement	62
3	For your safety	50	5.6	Transferring the measurement values via mobile network	62
3.1	Symbols used in this manual	50	6	Memory	63
3.2	Intended use	51	7	Error messages	65
3.3	Improper use	51	7.1	Blood pressure measurement errors	65
3.4	Children and persons with limited cognitive abilities	52	7.2	Communication errors	67
3.5	Batteries	52	8	Care and maintenance	69
3.6	Electrical fields	52	8.1	Care	69
3.7	Operation	53	8.2	Metrological monitoring	69
3.8	Ambient conditions	54	9	Disposal	70
3.9	Cleaning and disinfection	54	10	Technical data and symbols	71
4	Preparing the measurement	56	11	Standards	75
4.1	Unpacking	56	12	Warranty and repairs	75
4.2	Inserting the batteries	56	13	Manufacturer's EMC guidelines	77
			14	Frequency bands	82

Introduction

1 Introduction

Thank you for choosing the Tel-O-Graph® GSM or Tel-O-Graph® GSM plus upper arm blood pressure monitor.

Read these operating instructions carefully before use and keep them in a suitable place so that the information is available when required.

The Tel-O-Graph® GSM, GSM plus is a fully automatic blood pressure and pulse monitor. You have the option of transferring the data via a mobile network. The Tel-O-Graph® GSM plus also records the pulse wave form of the pulse and this information is transmitted alongside the blood pressure.

If you have any questions about services or products, feel free to contact us.

Contact data

Address	I.E.M. Industrielle Entwicklung Medizintechnik GmbH Cockerillstr. 69 52222 Stolberg Germany
Telephone	+49 2402 9500-93
Fax	+49 2402 9500-11
E-mail	service@iem.de
Website	www.iem.de

2 Description of device

2.1 The monitor


Fig. 1: The monitor

- 1 Display
- 2 Infrared interface (for service)
- 3  - button

- 4 Battery cover
- 5 Air tube socket

Description of device

2.2 Cuff


Fig. 2: Cuff

- 1 Cuff
- 2 Air tube
- 3 Air tube connection

2.3 Display


Fig. 3: Display

1	When this appears: Battery empty	6	Infrared
2	Display of systolic (upper) value	7	Number of pulse beats per minutes
3	Mobile network communication	8	Pulse detected
4	Signal strength	9	Number of measured values
5	Data transmission	10	Display of diastolic (lower) value

For your safety

2.4 Content

Content

- Blood pressure monitor
- Arm cuff (see label)
- Operating instructions
- Batteries (4x, AA, alkaline)

Accessories

- Arm cuff „S“ (arm circumference: 20-24 cm)
- Arm cuff „M“ (arm circumference: 24-32 cm)
- Arm cuff „L“ (arm circumference: 32-38 cm)
- Arm cuff „XL“ (arm circumference: 38-55 cm)

2.5 Clinical Validation

The accuracy of the device's measurements has been certified in accordance with ISO 81060-2:2013.

3 For your safety

This section explains all the safety aspects of the device.

Read this section carefully before using the monitor.

Contact your doctor before using the device if you are pregnant, are taking medication to alter blood clotting or if you have been diagnosed with cardiac arrhythmia or arteriosclerosis.

3.1 Symbols used in this manual

WARNING!

... indicates a potentially hazardous situation that could result in serious injury if not avoided.

i

... provides additional information or instructions on a section of this manual that should be considered.

3.2 Intended use

The Tel-O-Graph® GSM, GSM plus is used for the home measurement of the blood pressure and pulse on the upper arm.

The Tel-O-Graph® GSM plus additionally records data of the pulse wave.

The device is suitable for persons with an arm circumference of 20–55 cm when used with the corresponding cuff size.

It is designed for the transmission of the measured data via mobile network.

WARNING!

Self-diagnosis and self-treatment on the basis of the results is hazardous!

- Do not undertake any treatment and/or take medication due to the measured values without consulting your doctor.
- Follow your doctor's instructions.

Principal function

The principal function of the monitor is non-invasive measurement of blood pressure within the specified margin of error. If this is not possible, an error is displayed.

3.3 Improper use

The monitor must not be used for newborns, for operations, in the vicinity of MRI scanners or other strong magnetic fields, or for clinical monitoring of patients and during transport of patients.

i

The device has not been tested for use with pregnant women or those suffering from preeclampsia.

If you are taking medication to alter blood clotting, consult your doctor before using the monitor.

For your safety

3.4 Children and persons with limited cognitive abilities

Persons (including children) who are not able to use the monitor safely due to their physical, sensory or mental capabilities or their inexperience or lack of knowledge must not use this monitor without supervision or instruction by a responsible person. Keep device out of the reach of children and persons with limited cognitive abilities.

3.5 Batteries

- Remove the batteries from the battery compartment when they no longer have any charge or if you do not expect to use the monitor for a longer period of time.
- The batteries should not be exposed to fire or high temperatures.
- Do not attempt to recharge the batteries. Do not attempt to open or short-circuit the batteries. There is a risk of explosion.

- Liquid that escapes from the batteries due to mishandling can cause skin irritation. If you come into contact with the liquid, rinse it away with plenty of water. If the liquid comes into contact with your eyes, do not rub your eyes but instead immediately rinse them with water for 10 minutes and contact a doctor without delay.

3.6 Electrical fields

The measurements may be faulty if the device is operated in the vicinity of strong electrical fields. Do not operate the monitor near:

- High-voltage power lines
- Microwave devices

Portable and mobile HF transmitters, such as mobile phones, may affect the monitor.

Transmission of data via mobile communication networks may be disrupted by other devices, even if those devices correspond to the applicable transmission requirements specified by CISPR.

3.7 Operation

- Only use the cuffs approved for use with this monitor. Use of other cuffs or other parts not produced by the manufacturer may lead to incorrect measurement results.
- Do not drop the monitor and do not place objects on top of it.
- Do not place the cuff or the air tube around your neck when measuring. This can result in strangulation.
- Do not kink, knot or stretch the air tube for the cuff.
- Continuous pressure from the cuff can result in impaired blood circulation and injuries.
- Stop the measurement and remove the cuff if you feel any pain. Otherwise, this could lead to haematoma or bleeding.
- When a defibrillator is used, the monitor's housing must not touch the patient. This can result in damage and faulty readings.
- The monitor does not function for patients who suffer from atrial fibrillation.
- Do not place the cuff over wounds as this may result in injury.
- Do not place the cuff on an arm that has an intravascular catheter or arteriovenous fistula. This can result in temporary disruption of the blood flow.
- Do not take measurements too frequently. This can disrupt the blood flow.
- Consult your doctor if you wish to use the cuff on your arm on the side on which you have had a mastectomy.
- The monitor is subject to special precautionary measures with regard to EMC. For this reason, you must install and operate the device in accordance with the instructions specified in this manual.
- The monitor must not be placed directly next to or in a stack with other devices. If it is necessary to operate the device in a stack with other devices, examine the monitor to check that it is used as intended.

For your safety

3.8 Ambient conditions

- Never store the monitor outside a temperature range of -25 °C to +70 °C.
- Never store the monitor outside a temperature range of +5 °C to +40 °C.
- Only store or use the monitor at an air humidity (not condensing) of 15% to 93%.
- Extreme temperatures, humidity (not condensing) or altitudes can impact the performance of the monitor.

3.9 Cleaning and disinfection

- Only clean the outside of the monitor with a damp cotton cloth. Dampen the cotton cloth with lukewarm water only, possibly adding a mild detergent. Never use abrasive or solvent-based detergents as these could damage the surface of the monitor.
- Never use the monitor in the shower. If liquid has entered into the monitor, immediately remove the batteries from the monitor and send the device to your retailer or the manufacturer.
- The user (doctor) decides whether and when the cuff sleeve needs to be disinfected for hygienic reasons.

- I.E.M. GmbH recommends the following agents for disinfecting the cuff sleeve following successful inspection:
 1. Terralin Liquid (Manufacturer: Schülke & Mayr)
 2. Isopropanol (70%)

i

*Always observe the manufacturer's instructions on the use of these products.
Always allow the agents to dry without any residues.*

- Avoid disinfecting the cuff bladder and connected rubber tubing. Carefully remove the bladder (including the tubing) from the cuff sleeve before disinfection. The bladder and tubing can be damaged by disinfectants. Do not submerge the cuff sleeve in disinfectant solutions, as doing so may damage it.
- To achieve the optimum effect, apply the disinfectant to the cuff sleeve for at least five minutes.

- When using other disinfectants that are not recommended by I.E.M., the burden of proof that these can be used without causing damage lies with the user. Never use disinfectants which leave residues on the product or which are not suitable for contact with the skin.

i

There are patients who have intolerances (e.g. allergies) to disinfectants or its components. If you have such patients or are not sure if you have such patients, remove any potential residues by careful washing. Wash the cuff sleeve with a mild detergent in the washing machine at max. 30°C (without spin drying). Ensure that any used disinfectants are washed off, leaving no residues, before the blood pressure sleeve is worn.

- Do not use fabric softeners or other additives (e.g. hygiene rinses, textile deodorants). These agents can leave behind residues and damage the material. The cuff sleeve is not suitable for drying in a dryer. The Velcro strip must always be closed before washing.

Preparing the measurement

4 Preparing the measurement

4.1 Unpacking

WARNING!

Risk of strangulation from the air tube and cuff!

- Keep the cuff out of the reach of children!

- 1) Unpack everything delivered and check that everything is there.
- 2) Inspect the monitor to ensure there is no visible damage on the outside. If the monitor is damaged, have it repaired before use.
- 3) Keep the packaging so that it can be safely packaged at a later date.

4.2 Inserting the batteries


Fig. 4: Opening the battery cover

- 1) Open the battery cover on the underside of the monitor.
- 2) Insert four batteries (AA, alkaline), with the poles aligned correctly.
- 3) Close the battery cover.

⇒ The monitor switches on automatically.

4.3 Switching the monitor on/off

Switching on

Fig. 5: Start-screen

- 1) Press the  - button.

⇒ The display shows the start screen (Fig. 5).

Switching off

- 2) The monitor switches off automatically after 5 minutes.

Measuring blood pressure and pulse

5 Measuring blood pressure and pulse

5.1 Before measuring

- Choose the right cuff size.
- Avoid eating, smoking or any strenuous activity directly before the measurement. All these factors affect the results. Try to relax in a calm atmosphere for about 5 minutes before measuring your blood pressure.
- Always measure on the same arm (normally your left).
- Try to carry out the measurements at regular intervals and at the same time of day as your blood pressure changes over the course of the day.

i

Larger or smaller cuffs are available via the service address.

5.2 Armmanschette anlegen

WARNING!

Risk of injury due to incorrectly connected air tube!

- Do not kink, knot or stretch the air tube.

WARNING!

The monitor must only be operated with the original cuff as there is otherwise a risk of incorrect readings or injury.


Fig. 6: Insert air tube

- 1) Insert the air tube connection into the air tube socket on the left-hand side of the monitor (Fig. 6).
- 2) Uncover your upper arm.

i

The cuff must be directly against your skin.

- 3) Insert your left arm into the cuff.


Fig. 7: Correct position of cuff

- 4) Position the cuff correctly (Fig. 7):
 - The tube must run down the middle of the inside of your lower arm towards your middle finger.
 - The lower edge of the cuff must be about 2 cm above your elbow.
 - The artery marking (arrow above the tube) must be pointing downwards on the middle of the inside of your arm.
- 5) Once the cuff is positioned correctly, fasten the end with the Velcro.

Measuring blood pressure and pulse

i

Place the cuff so that there is enough room to insert your index and middle fingers between the skin and the cuff's material.

5.3 Correct posture

After positioning the cuff, assume the correct posture to take the measurements (Fig. 8).


Fig. 8: Correct posture

- 1) Sit down comfortably on a chair to measure your blood pressure.
- 2) Place your elbow on a table or a solid surface.

- 3) Hold your arm so that the cuff is at the same height as your heart.
- 4) Relax your arm and turn the palm of your hand to face upwards.
- 5) Place your feet flat on the floor and leave your legs uncrossed.

i

The reading can be affected by the position of the cuff, the patient's correct posture and position, physical exertion or the patient's mental state.

5.4 Measuring

WARNING!

Risk of injuries during measurement!

- Do not place the air tube around your neck.
- If you feel any pain during the measurement, stop the measurement.
- Do not place the cuff over wounds.

- 1) Apply the cuff and assume the correct posture.
- 2) Press the  – button.
⇒ The display shows the start screen (Fig. 5).
- 3) Press the  – button again to start the measurement.
⇒ The monitor confirms this with a short acoustic signal. The cuff slowly pumps up. After the initial pumping, there is additional pumping. Once the monitor detects a pulse, the icon ♥ appears. The reading is taken as the air is let out.

i

During a pulse reading with the Tel-O-Graph® GSM plus, the cuff pressure will be kept constant during the air release for 10 seconds at the diastolic blood pressure value to enable the pulse wave form to be read.

i

Do not talk during the measurement.

You can stop the measurement at any time by pressing the  – button.

After the measurement, all the air is let out of the cuff. The display shows your blood pressure and your pulse rate.

- 4) Check the results.
- 5) Remove the cuff.

Measuring blood pressure and pulse

WARNING!

Taking measurements too frequently can result in injury!

Wait for at least 2 or 3 minutes to give your arteries time to return to their initial state.

5.5 Stopping the measurement

- 1) If the measurement causes you pain, press the

– button during the measurement.

⇒ The air is let out of the cuff and the measurement is cancelled.

- 2) Remove the cuff.

5.6 Transferring the measurement values via mobile network

If the monitor is set up for transmission of your measurement values via a mobile network, the measurement values will be automatically transferred to the database.

If previous measurement transmissions have not yet been transferred, all data not yet sent will be sent together.

The monitor emits a short acoustic signal to confirm successful transmission of the reading.

6 Memory

Saving readings


Fig. 9: Number of readings

The monitor can save 15 blood pressure readings and 15 pulse rate readings.

i

Only readings that have not yet been transferred to the database are saved.

If there are more than 15 readings, the oldest data is overwritten by the new readings.

i

The data in the database is retained and is not overwritten.

The number of readings is shown on the display (Fig. 9).

Deleting readings from the monitor

To clear the memory on your monitor, proceed as follows:


Fig. 10: Delete measurements

- 1) Press and hold the  – button on the monitor for about 5 s.

⇒ The display shows [IP] (Fig. 10/A).

i

The display automatically shows each menu entry for 5 seconds. The display order is: Mobile network transmission [IP] → Infrared transmission [IR] → Delete measurement values (c Ir)

Memory

2) Wait until *[c lr]* flashes on the display (Fig. 10/B).

3) Press the  – button.

⇒ You hear an acoustic signal and the display shows *[c lr]* permanently.

4) Press and hold the  – button for more than 5 s.

⇒ The acoustic signal sounds 3 times and the display shows *[M00]* (Fig. 10/C).

The measurements are deleted.

7 Error messages

7.1 Blood pressure measurement errors

The monitor indicates blood pressure measurement errors and communication errors by emitting 12 short acoustic signals.

EN

Fault description	Cause	Remedy
Err 1	Arm moved during measurement.	Keep your arm still during the measurement (→ Chapter 5.3).
	Insufficient valid pulse rate detected.	Reposition the cuff (→ Chapter 5.2).
Err 2	Arm moved during measurement.	Keep your arm still during the measurement (→ Chapter 5.3).
	The cuff is not correctly positioned.	Check that the cuff is fitted properly (→ Chapter 5.2).
Err 3	Blood pressure beyond the measurement range.	If this notification appears continuously (or repeatedly), it is possible that the monitor is not suitable for you. Contact your doctor.
	Strong arm movement.	Keep your arm still during the measurement (→ Chapter 5.3).

Error messages

Fault description	Cause	Remedy
Err 4	Readings not transmitted correctly.	Wait until you take another measurement. The data are then transmitted together with the new readings. If this error occurs repeatedly, send the monitor to your retailer or directly to the manufacturer for checking.
Err 5 	Power pack or battery voltage too low.	Change the batteries (→ Chapter 4.2).
	Battery contacts corroded.	Clean the battery contacts with a cotton cloth and a little alcohol.
Err 6	Air congestion.	Check the cuff for air congestion or a kink in the air tube. If there is a kink in the air tube, unkink it otherwise send the monitor to your retailer or the manufacturer.
	Cuff not connected properly.	Connect the cuff to the monitor (→ Chapter 5.2).
	Leak in the cuff or air tube.	Replace the cuff and air tube.
Abr. (Abort)	Measurement cancelled on request.	
Err 9	Residual pressure in the cuff.	Wait until the cuff has vented completely.

Fault description	Cause	Remedy
Err 10	Internal error.	If this error occurs repeatedly, send the monitor to your retailer or directly to the manufacturer for checking.

7.2 Communication errors

Fault description	Cause	Remedy
Cod 1	Registration with mobile network not possible.	Move to a location with better mobile network reception. If this error occurs repeatedly, contact the specialist dealership.
Cod 2	No GSM connection.	Move to a location with better mobile network reception. If this error occurs repeatedly, contact the specialist dealership.
Cod 3	Could not send readings.	Take another measurement and check whether the readings are transmitted. If this error occurs repeatedly, contact the specialist dealership.
Cod 4	Communication error.	If this error occurs repeatedly, contact the specialist dealership.
Cod 5	Invalid target data in EEPROM (memory).	If this error occurs repeatedly, contact the specialist dealership.

Error messages

Fault description	Cause	Remedy
Cod 6	Hardware communication error.	If this error occurs repeatedly, contact the specialist dealership.
Cod 7	Memory empty.	If this error occurs repeatedly, contact the specialist dealership.
Cod 8	No PWA license available.	If this error occurs repeatedly, contact the specialist dealership.
Cod 9	SIM card not inserted.	Insert the SIM card for the device.
Cod 10	Other error.	If this error occurs repeatedly, contact the specialist dealership.

8 Care and maintenance

8.1 Care

**WARNING!**

Risk of damage from water!

- Switch the device off and remove the batteries.

If necessary, wipe the outside only of the monitor and cuff with a damp cotton cloth. Dampen the cotton cloth with lukewarm water only, possibly adding a mild detergent. Never use chemical solvents and detergents as these could damage the surface of the monitor.

If, despite this, water enters the monitor, contact your retailer or the manufacturer without delay.

Remove the batteries if you are not going to use the monitor for a longer period of time.

8.2 Metrological monitoring

The monitor and cuff are calibrated by the manufacturer for a period of two years. The metrological check must be carried out in accordance with directive 93/42/EEC every two years at the latest if the device is used professionally. In certain countries, this requirement may be regulated by national laws or regulations.

The metrological check is subject to a charge and can be carried out by **I.E.M. GmbH**, by a responsible public authority or by an authorized maintenance service in accordance with the "*Medical Devices Operator Ordinance*".

Disposal

9 Disposal

Monitor

The symbol on the product or packaging means that this product should be treated as normal domestic waste, but has to be passed to a recycling point for electric and electronic devices.

You can find out more about this from your local authority, the communal disposal companies or the shop in which you bought the product.

Batteries

Batteries must not be disposed of with domestic waste. As a consumer you are legally bound to give back used batteries. You can return your old batteries to official collection points within your community or everywhere batteries of that type are sold.

	Li	Battery contains lithium
	Al	Battery contains alkali
	Mn	Battery contains manganese

10 Technical data and symbols

Technical data

The monitor complies with the EMC directives.

The cuff and air tube are made of a non-conductive material. This makes the device defibrillator-safe.

EN

Specification	Value	Unit
Measuring method	Oscillometric	
Blood pressure measurement range	30 to 290	mmHg
Pulse measuring range	30 to 240	1/Min
Pressure accuracy	+/- 3	mmHg
Memory	15	Measurements
Power supply	6V (4x AA, 1,5V, alkaline)	
Battery capacity	200	Measurements
Data connection	Mobile network	
Data transmission	GPRS Class 12	
Dimesions (L x W x H)	150 x 109 x 56	Mm
Weight (without batteries)	325	g

Technical data and symbols

Specification	Value	Unit
Material	ABS (Acrylnitril-Butadien-Styrol)	
IP protection class	21	
Operating temperature	+5 to +40	°C
Ambient pressure	700 to 1060	hPa
Transport temperature	-25 to +70	°C
Storage temperature	-25 to +70	°C
Humidity, not condensing (operation, transport and storage)	15 to 93	%

Symbols

Symbol	Meaning
	4 x LR6 or HR6, AA with 1.5 V
	Manufacturer in accordance with Directive 93/42/EEC
	Designation of the FCC radio authorization
	CE 0044: Designation of a medical device in accordance with Directive 93/42/EEC
	It may not be possible to use the device everywhere in Europe.
	Protect from rain and moisture
	Comply with the operating instructions
	The symbol on the product or packaging means that this product should not be treated as normal domestic waste, but should be passed to a recycling point for electric and electronic devices. You can find out more about this from your local authority, the communal disposal companies or the shop in which you bought the product.

Technical data and symbols

Symbol	Meaning
	Device is defibrillator-safe
	Non-ionising radiation

11 Standards

The device complies with the specifications of EU Directive 93/42/EEC for medical products and with the national laws on medical products.

The device complies with the requirements of Directive 1999/5/EC (R&TTE). The monitor visibly bears the CE mark.

The monitor fulfills Directive 2002/95/EC of the European Parliament and the European Council of 27 January 2003 on restricting the use of specific hazardous materials in electrical and electronic equipment.

12 Warranty and repairs

Warranty information

- **I.E.M. GmbH** provides a two year warranty on the blood pressure monitor from the date of sale. Proof of the date of sale shall be provided in the form of a properly completed warranty card or an invoice.
- Faults due to material or production defects shall be remedied free of charge within the warranty period.
- A warranty claim does not result in an extension of the warranty period, neither for the device nor for the replaced components.
- The following are excluded from the warranty:
 - All damage caused by improper handling, e.g. by non-compliance with the operating instructions.
 - Damage resulting from maintenance or intervention by the purchaser or unauthorized third parties.
 - Transport damage on route from the manufacturer to the consumer or during shipment to the customer service.
 - Accessory parts subject to normal wear and tear (cuff, batteries, etc.).

Warranty and repairs

- Liability for direct or indirect consequential damages caused by the device is excluded even in the event that the damage to the device is acknowledged as a warranty case.
- Warranty cards that have been amended without authorization are invalid.
- Any further claims, irrespective of the cause, are excluded.

I.E.M. GmbH does not grant any warranty on the batteries provided.

Repair conditions

In the event of a warranty claim, fill out the warranty card with the precise information about the sender and fault. Please contact our customer service, who will inform you of the shipping arrangements.

13 Manufacturer's EMC guidelines

Electromagnetic emissions

The Tel-O-Graph® GSM, Tel-O-Graph® GSM plus is designed for use in the environment specified below. Only use the Tel-O-Graph® GSM, Tel-O-Graph® GSM plus in such an environment.

EN

Emitted interference measurement	Compliance	Electromagnetic environment - guidelines
HF Emissions according to CISPR 11	Group 1	The Tel-O-Graph® GSM, Tel-O-Graph® GSM plus only uses HF energy for its internal functionality. Its HF emission is therefore very low and it is improbable that any neighboring electronic device will experience any interference.
HF Emissions according to CISPR 11	Class B	The Tel-O-Graph® GSM, Tel-O-Graph® GSM plus is suitable for use in other facilities than the living area and those immediately connected to the public supply network, which also supplies buildings used for residential purposes.
IEC 61000-3-2	Not applicable	
IEC 61000-3-3	Not applicable	

Manufacturer's EMC guidelines

Electromagnetic immunity

The Tel-O-Graph® GSM, Tel-O-Graph® GSM plus is designed for use in the environment specified below. Only use the Tel-O-Graph® GSM, Tel-O-Graph® GSM plus in such an environment.

Immunity test	IEC 60601-1-test level	Compliance level	Electromagnetic environment - guidelines
Electrostatic discharge (ESD) in accordance with IEC 61000-4-2	+/- 6kV contact discharge +/- 8kV air discharge	+/- 6kV contact discharge +/- 8kV air discharge	Floors should consist of wood or concrete or be tiled with ceramic tiles. If the floor consists of synthetic materials, relative humidity must be at least 30 %.
IEC 61000-4-4		Not applicable	
IEC 61000-4-5		Not applicable	
IEC 61000-4-11		Not applicable	
Magnetic field at supply frequency (50/60 Hz) in accordance with IEC 61000-4-8	3 A/m	3 A/m	Magnet fields in network frequency should correspond with the typical values found in business and hospital environments.

Immunity test	IEC 60601-1-test level	Compliance level	Electromagnetic environment - guidelines
Radiated RF disturbance in accordance with IEC 61000-4-3	3 V/m 80 Mhz to 2,5 Ghz	3 V/m	Portable and mobile wireless receiving and transmission equipment should not be used at a distance to the Tel-O-Graph® GSM, Tel-O-Graph® GSM plus (including wires) of less than the recommended protection ratio calculated for the equation that applies to the transmission frequency. Recommended safety distance: $d = 1,2 * \sqrt{P(W)}$ for 80 Mhz to 800 Mhz $d = 2,3 * \sqrt{P(W)}$ for 800 Mhz to 2,5 Ghz where P = maximum nominal output of the transmitter in watts as per the transmitter manufacturer's specifications and d = recommended protection ratio in meters (m). The field strength of stationary radio wireless transmitters should be lower than the compliance level at all frequencies, as shown by an on-site examination. Interference is possible in the vicinity of devices bearing the graphic symbol "nonionising radiation".

Manufacturer's EMC guidelines

NOTE 1: The higher frequency range applies for 80 MHz and 800 MHz.

NOTE 2: These guidelines may not apply in all cases. The diffusion of electromagnetic parameters is affected by absorption and reflections from buildings, objects and people.

The field strength of stationary transmitters, such as base stations for cordless telephones and land mobile radios, amateur radio stations, AM and FM radio and television transmitters cannot be specified precisely in advance, in theory.

In order to determine the electromagnetic environment relating to stationary transmitters, you should consider a study of the location. If the measured field strength at the location where the Tel-O-Graph[®] GSM, Tel-O-Graph[®] GSM plus is used exceeds the upper compliance level, the Tel-O-Graph[®] GSM, Tel-O-Graph[®] GSM plus should be kept under observation to ensure proper functioning. If unusual performance characteristics are observed, additional action may be necessary, such as changing the orientation or location of the Tel-O-Graph[®] GSM, Tel-O-Graph[®] GSM plus.

Recommended safety distances between portable and mobile HF telecommunication devices and the Tel-O-Graph® GSM, Tel-O-Graph® GSM plus

The Tel-O-Graph® GSM, Tel-O-Graph® GSM plus is intended for operation in an electromagnetic environment in which RF disturbance is under control. Keep the minimum distances between portable and mobile RF telecommunication devices (transmitters) and the monitor specified below, depending on the output of the communication device.

EN

Protection ratio in relation to transmission frequency in m			
Nominal output of transmitter P in W	80 Mhz to 800 Mhz $d = 1,2 * \sqrt{P(W)}$	800 Mhz to 2,5 Ghz $d = 2,3 * \sqrt{P(W)}$	150 Khz to 80 Mhz
0,01	0,12	0,23	Not applicable
0,1	0,38	0,73	Not applicable
1	1,2	2,3	Not applicable
10	3,8	7,3	Not applicable
100	12	23	Not applicable

The recommended protection ratio can be calculated in meters (m) for transmitters whose maximum nominal output is not specified in the table above. To do so, use the equation for the column in question, where P is the maximum nominal output of the transmitter.

NOTE 1: The higher frequency range applies for 80 MHz and 800 MHz.

NOTE 2: These guidelines may not apply in all cases. The diffusion of electromagnetic parameters is affected by absorption and reflections from buildings, objects and people

Frequency bands

14 Frequency bands

The Tel-O-Graph[®] GSM, Tel-O-Graph[®] GSM plus uses the following frequency bands:

Band designation	Range	Transmission class
E-GSM	880 MHz to 915 MHz	Class 4 (+33dBm $\pm$ 2dB)
DCS	1710 MHz to 1785 MHz	Class 1 (+30dBm $\pm$ 2dB)
FDD I	1920 MHz to 1980 MHz	Class 3 (+24dBm +1/-3dB)
FDD VIII	880 MHz to 915 MHz	Class 3 (+24dBm +1/-3dB)

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La société **I.E.M. GmbH** n'accorde aucune garantie sur les piles fournies.